



COMMERCIAL RESPIRATOR 95%

INTENDED USE

The Commercial Respirator 95% is intended to be worn to protect both the patient and healthcare personnel from the transfer of microorganisms, body fluids and particulate material. These face masks are intended for use in infection control practices to reduce the potential exposure to blood and bodily fluids. This is a single use, disposable device.



PRODUCT INFORMATION

SPECIFICATIONS:

- Testing per 42CFR84
- NIOSH rating
- Filters 95% or more of relevant particles
- PFE, Particulate Efficiency >95% at 0.0725 microns
- BFE, Bacteria Live Particle Filtration Efficiency >95% 0.1 micron
- Flame Spread Class 1
- ASTM F2100-11 Standard ASTM Level 1

MATERIALS: 4 plies of fabric consisting of:

- 2 layers spun bonded polypropylene
- 2 layers melt blown nonwoven polypropylene
- Twist tie wire sonic welded inside the mask
- Elastic straps to be placed around the head

1. Product is valid for 2 years. Please use within 2 years of manufacture date.
2. This product is for single-use only. Please discard after each use.

INSTRUCTIONS FOR USE

1. Open package and remove mask. When wearing the mask, ensure the white side is facing inwards, and metal nose clip is facing up.
2. Place mask against face and wrap loops around ears. Gently pinch the metal nose clip against the bridge of the nose to secure in place.

NOTE: This product is for one-time use. Dispose of it according to the requirements of medical waste management.

Authorized Distributor:
DISK-O-TAPE, INC.

phone | 216-765-8273
www.disk-o-tape.com

23775 Mercantile Road
Cleveland, Ohio 44122-5917